

## The more you know...

In TogetherCare, regular insulin sliding scales for our patients receiving enteral feeds or TPNs, can only be found in the **enteral feed** order set.


1. To find the closest option to our current orders, select the Basal and Correctional insulin options highlighted



2. Next, order your long acting insulin, then select "correctional insulin regular"



3. Based on your patients TDD insulin requirements, select the appropriate scale



Please note, there is no variation for the rate of tube feeds in TC.

They are also considerably less aggressive than our current enteral feed scales.

## Our current SJH scales, for reference

## **Continuous Enteral Feeding**

Supplemental Insulin Scale – **Regular Insulin** SQ every 6 hours while enteral feeding is running

Category	□ Frail		□ Average		□ Obese		Custom
Enteral Feed Rate	<50 (ml/hr)	≥ 50 (ml/hr)	<50 (ml/hr)	≥ 50 (ml/hr)	<50 (ml/hr)	≥ 50 (ml/hr)	Scale
Blood glucose (mg/dl)							
121-170	2 units	4 units	5 units	7 units	10 units	12 units	units
171-220	3 units	5 units	7 units	9 units	12 units	14 units	units
221-270	4 units	6 units	8 units	10 units	14 units	16 units	units
271-320	5 units	7 units	10 units	12 units	18 units	20 units	units
321-370	6 units	8 units	13 units	16 units	20 units	24 units	units
371-420	6 units	8 units	14 units	17 units	22 units	26 units	units
>420	7 units	9 units	15 units	18 units	24 units	28 units	units