

St. Joseph's Health Syracuse Workflow: Lab Specimen Labels: Which Labels to Order in PeopleSoft



Upon the Wave 8 Epic Beaker Lab go-live, each office and IP unit will be responsible for ordering and maintaining stock of specimen label printer label stock to be used in the printers used in the department. Below is a reference detailing the PeopleSoft products for the two primary specimen label printers that have been deployed.

NOTE: Please ensure adequate supplies are available based on your anticipated specimen volume for your department. The Laboratory will not be able to provide these labels if you run out.

For Zebra ZQ610 Printers (Grab & Go Printers, Mobile) – Primarily Acute/IP Areas



These printers are the smaller, mobile printers used in Med-Surg areas as “grab-and-go” printers. They are also located in Critical Care and Emergency Department and associated with specific WOW’s.

MFR# TD75F-1142 PS# 726781 Direct Thermal Synthetic WHITE

For Zebra ZQ621 Printers (Desktop printers) – Primarily in Ambulatory and OR



These printers are larger desktop printers. They are primarily located in Ambulatory Offices and the OR and procedural areas.

MFR# LDWHP173 PS#363814 Synthetic

***For questions related to PeopleSoft, please contact your Supply Chain Management resource.